



Socks Order Form

Order No.: _____ Patient Reference No.: _____

First select the required garment, then tick the box to confirm the standard configuration, or choose from the style options on page 4:

Garment (please indicate)

☐ PO 0020 Anklet ☐ Left ☐ Right 	☐ PO 0021 Knee Length Sock ☐ Left ☐ Right 
---	--

Standard configurations for Socks



Feature	☐ Anklet & Knee Length Sock (0020,0021)
Proximal	Elastic cuff 2.5cm
Distal	Overlock / open toes

Please note: you may use one form to order both left and right socks, however they must be the same garment type and theme/fabric

Socks Order Form

Order No.: _____ Patient Reference No.: _____

Please use this outline in conjunction with the guide and grids below
Use a traditional tape measure to record measurements. Type them in the boxes below or in the corresponding grid on the next page, both fill simultaneously.



Circumference Measurements		Left (cm)	Right (cm)
15	Knee joint		
16	Upper calf		
17	Mid calf		
18	Lower calf		
19	Upper margin of medial malleolus		
24	Around foot and heel under malleolus		
25	Instep or waist of foot		
26	Foot at metatarsal heads		

Length Measurements		Left (cm)	Right (cm)
21	Medial malleolus to required top of sock		
27	Upper margin medial malleolus to sole of foot		
28	Metatarsal heads to heel (medial)		
28a	From metatarsal heads to heel on lateral border		
29	Metatarsal heads to tip of toes (medial)		

Style options when standard configuration is not selected		Left	Right
Proximal elastic:	Overlock (no elastic)	☐	☐
	Regular (inverted) 1.5cm	☐	☐
	Regular (inverted) 2.5cm	☐	☐
	Cuff 1.5cm	☐	☐
	Cuff 2.5cm	☐	☐
	Silicone Regular (inverted) 2.5cm	☐	☐
	Silicone cuff 2.5cm	☐	☐
Ankle contracture seam (at front of ankle for shaping only)		☐	☐
Distal elastic:	Overlock (no elastic)	☐	☐
	Regular (inverted) 1.5cm	☐	☐
	Regular (inverted) 2.5cm	☐	☐
	Cuff 1.5cm	☐	☐
	Cuff 2.5cm	☐	☐
	Silicone Regular (inverted) 2.5cm	☐	☐
	Silicone cuff 2.5cm	☐	☐

Socks Order Form

Order No.: _____ Patient Reference No.: _____

Modifications available for all garments

All the following items will be an additional charge

Zippers - 1145 (tick if required)

Zipper placement	☐ Inside of fabric	☐ Outside of fabric
Position (please select: Medial, Lateral or Posterior)	Left	Right
	☐	☐
Length: _____ cm		
Hook & eye (on fly behind the zip to assist donning)	☐	☐

Inset Zippers - 1144 (tick if required)

Use placement pad to mark position if required		
Zipper placement	☐ Inside of fabric	☐ Outside of fabric
Position (please select: Medial, Lateral or Posterior)	Left	Right
	☐	☐
Length: _____ cm		

Reinforcements

Reinforced heel (for high wear area to reinforce)	1187	☐	☐
Non-slip silicone sole of foot	1188	☐	☐

Toes

Self enclosed toe in base fabric (no seams)	1159	☐	☐
Soft enclosed toe in lining fabric	1160	☐	☐

Silon-TEX® II Insert

Silon-TEX® II fabric (sewn into garment)	1191	☐ Use placement pad to mark position
--	------	---

Pockets & Pads

Pocket (sewn in for the insertion of pads to apply extra pressure to certain areas) Please specify position	0027	☐ Use placement pad to mark position
Silon-TEX® II pocket (as above but covered with Silon-TEX® II fabric)	1147	☐ Use placement pad to mark position

Foam Pads (to insert into pocket, please select foam thickness)

Low profile 5mm	1178	☐
Low density 20mm	1179	☐